



Recombinant Human Zinc finger protein 670 (ZNF670)

Product Code	CSB-MP026936HU
Abbreviation	ZNF670
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BS34
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSVSFEDVA VAFTQEEWAL LDPSQKNLYR DVMQEIFRNL ASVGNKSEDQ NIQDDFKNPG RNLSSHVVER LFEIKEGSQY GETFSQDSNL NLNKKVSTGV KPCECSVCGK VFICHSALHR HILSHIGNKL FECEECPEKL YHCKQCGKAF ISLTSVDRHM VTHTSNGPYK GPVYEKPFDF PSVFQMPQST YTGEKTYKCK HCDKAFNYSS YLREHERTHT GEKPYACKKC GKSFTFSSSL RQHERSHTGE KPYECKECGK AFSRSTYLG I HERTHTGEKP YECIKCGKAF RCSRVLRVHE RTHSGEKPYE CKQCGKAFKY SSSLCEHERT HTGVKPYGCK ECGKSFTSSS ALRSHERTHT GEKPYECKKC GKAFSCSSSL RKHERAYMW
Source	Mammalian cell
Target Names	ZNF670
Protein Names	Recommended name: Zinc finger protein 670
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.