



Recombinant Human RING finger protein 208 (RNF208)

Product Code	CSB-MP019871HU
Abbreviation	RNF208
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0X6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPSDPGPEAG SGWPGLLMSC LKGPHVILKM EAMKIVHPEK FPELPAAPCF PPAPRPTPTL APKRAWPSDT EIVNQACGG DMPALEGAPH TPPLPRRPRK GSSELGFPRV APEDEVIVNQ YVIRPGPSAS AASSAAAGEP LECPTCGHSY NVTQRRPRVL SCLHSVCEQC LQILYESCPK YKFISCPTCR RETVLFTDYG LAALAVNTSI LSRLPPEALT APSGGQWGAE PEGSCYQTFR QYCGAACTCH VRNPLSACSI M
Source	Mammalian cell
Target Names	RNF208
Protein Names	Recommended name: RING finger protein 208
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.