



Recombinant Human FGFR1 oncogene partner 2 (FGFR1OP2)

Product Code	CSB-MP008644HU
Abbreviation	FGFR1OP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NVK5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSCTIEKALA DAKALVERLR DHDDAAESLI EQTTALNKR V EAMKQYQEEI QELNEVARHR PRSTLVMGIQ QENRQIRELQ QENKELRTSL EEHQSALELI MSKYREQMFR LLMASKKDDP GIIMKLKEQH SKIDMVHRNK SEGFFLDASR HILEAPQHGL ERRHLEANQN ELQAHVDQIT EMAAVMRKAI EIDEQQGCKE QERIFQLEQE NKGLREILQI TRESFLNLRK DDASESTSL S ALVTNSDSL RKS
Source	Mammalian cell
Target Names	FGFR1OP2
Protein Names	Recommended name: FGFR1 oncogene partner 2
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.