



# Recombinant Human Interleukin-1 receptor-associated kinase 4 (IRAK4)

<b>Product Code</b>	CSB-EP011812HU-B
<b>Abbreviation</b>	IRAK4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NWZ3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNKPITPSTY VRCLNVGLIR KLSDFIDPQE GWKKLAVAIK KPSGDDRYNQ FHIRRFALL QTGKSPTSEL LFDWGTTNCT VGDLVDLLIQ NEFFAPASLL LPDAVPKTAN TLPSKEAITV QKQMPFCDK DRTLMTVPVQN LEQSYMPPDS SSPENKSLEV SDTRFHSFSF YELKNVTNNF DERPISVGGN KMGEGGFGVV YKGYVNNTTV AVKKLAAMVD ITTEELKQF DQEIKVMAC QHENLVELLG FSSDGDDLCL VYVYMPNGSL LDRLSCLDGT PPLSWHMRCK IAQGAANGIN FLHENHHIHR DIKSANILLD EAFTAKISDF GLARASEKFA QTVMTSRIVG TTAYMAPEAL RGEITPKSDI YSFGVVLLEI ITGLPAVDEH REPQLLLDIK EEIEDEEKTI EDYIDKKMND ADSTSVEAMY SVASQCLHEK KNKRPDIKKV QQLLQEMTAS
<b>Source</b>	E.coli
<b>Target Names</b>	IRAK4
<b>Protein Names</b>	Recommended name: Interleukin-1 receptor-associated kinase 4 Short name= IRAK-4 EC= 2.7.11.1 Alternative name(s): Renal carcinoma antigen NY-REN-64
<b>Expression Region</b>	1-460
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.