



# Recombinant Human Zinc finger protein 167 (ZNF167)

<b>Product Code</b>	CSB-MP026552HU
<b>Abbreviation</b>	ZNF167
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9P0L1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTTAGRGNLG LIPRSTAFQK QEGRLTVKQE PANQTWGQGS SLQKNYPPVC EIFRLHFRQL CYHEMSGPQE ALSRLRELCR WWLMPEVHTK EQILELLVLE QFLSILPGEL RTWVQLHHP E SGEEAVAVVE DFQRHLSGSE EVSAPAQKQE MHFEETTALG TTKESPPTSP LSGGSAPGAH LEPPYDPGTH HLP SGDFACQ TSPVPTLPQV GNSGDQAGAT VLRMVRPQDT VAYEDLSVDY TQKKWKS LTL SQRALQWNMM PENHHSMASL AGENMMKGSE LTPKQEFFKG SESSNRTSGG LFGVVP GAAE TGDVCEDTFK ELEGQTSDEE GSRL ENDFLE ITDEDKKST KDRYDKYKEV GEHPPLSSSP VEHEGVLKGQ KSYRCDECGK AFNRSSHLIG HQR IHTGEKP YECNECGKTF RQTSQLIVHL RTHTGEKPYE CSECGKAYRH SSHLIQHQL HNGEKPYKCN ECAKAFQSS RLTDHQRTH GEKPYECNEC GEAFIRSKSL ARHQVLHTGK KPYKCN ECGR AFCSNRNLID HQR IHTGEKP YECSECGKAF SRSKCLIRHQ SLHTGEKPYK CSECGKAFNQ NSQLIEHERI HTGEKPF ECS ECGKAFGLSK CLIRHQRLHT GEKPYKCN E GKSFNQNSHL I IHQR IHTGE KPYECNECGK VFSYSSSLMV HQRTH TGEKP YKCNDCGKAF SDSSQLIVHQ RVHTGEKPYE CSECGKAFSQ RSTFNHHQRT HTGEKSSGLA WSVS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZKSCAN7
<b>Protein Names</b>	Recommended name: Zinc finger protein 167 Alternative name(s): Zinc finger protein 448 Zinc finger protein 64 Zinc finger protein with KRAB and SCAN domains 7
<b>Expression Region</b>	1-754
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.