



# Recombinant Human Myotilin (MYOT)

<b>Product Code</b>	CSB-YP015364HU
<b>Abbreviation</b>	MYOT
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UBF9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFNYERPKEHF IQSQNPCGSR LQPPGPETSS FSSQTKQSSI IIQPRQCTEQ RFSASSTLSS HITMSSSAFP ASPKQHAGSN PGQRVTTTTYN QSPASFLSSI LPSQPDYNSS KIPSAMDNY QQSSAGQPIN AKPSQTANAK PIPRTPDHEI QGSKEALIQD LERKLCCKDT LLHNGNQRLT YEEKMARRLL GPQNAAAVFQ AQDDSGAQDS QQHNSEHARL QVPTSQVRSR STSRGDVNDQ DAIQEKFYPP RFIQVPENMS IDEGRFCRMD FKVSGLPAPD VSWYLNVRTV QSDDLHKMIV SEKGLHSLIF EVVRASDAGA YACVAKNRAG EATFTVQLDV LAKEHKRAPM FIYKPQSKKV LEGDSVKLEC QISAIPPPKL FWKRNNEMVQ FNTDRISLYQ DNTGRVTLI KDVNKKDAGW YTVSAVNEAG VTTCNTRLVD TARNQTLPA PKQLRVRPTF SKYLALNGKG LNVKQAFNPE GEFQRLAAQS GLYESEEL
<b>Source</b>	Yeast
<b>Target Names</b>	MYOT
<b>Protein Names</b>	Recommended name: Myotilin Alternative name(s): 57 kDa cytoskeletal protein Myofibrillar titin-like Ig domains protein Titin immunoglobulin domain protein
<b>Expression Region</b>	1-498
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.