



Recombinant Human Prenylcysteine oxidase 1 (PCYOX1)

Product Code	CSB-BP017652HU
Abbreviation	PCYOX1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UHG3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AEL RAPPDKIAII GAGIGGTSAA YYLRQKFGKD VKIDLFEREE VGGRLATMMV QGQEYEAGGS VIHPLNLHMK RFVKDLGLSA VQASGGLLGI YNGETLVFEE SNWFIINVIK LVWRYGFQSL RMHMMWVEDVL DKFMRIYRYQ SHDYAFSSVE KLLHALGGDD FLGMLNRTLL ETLQKAGFSE KFLNEMIAPV MRVNYGQSTD INAFVGAVSL SCSDSGLWAV EGGNKLVCSSG LLQASKSNLI SGSVMYIEEK TKTKYTGNTPT KMYEVVYQIG TETRSDFYDI VLVATPLNRK MSNITFLNFD PPIEEFHQYY QHIVTTLVKG ELNTSIFSSR PIDKFGLNTV LTTDNSDLFI NSIGIVPSVR EKEDPEPSTD GTYVWKIFSQ ETLTKAQILK LFLSYDYAVK KPWLAYPHYK PPEKCPSSIL HDRLYYLNIGI ECAASAMEMS AIAAHNAALL AYHRWNGHTD MIDQDGLYEK LKTEL
Source	Baculovirus
Target Names	PCYOX1
Protein Names	Recommended name: Prenylcysteine oxidase 1 EC= 1.8.3.5 Alternative name(s): Prenylcysteine lyase
Expression Region	28-505
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.