



# Recombinant Human Prenylcysteine oxidase 1 (PCYOX1)

<b>Product Code</b>	CSB-MP017652HU
<b>Abbreviation</b>	PCYOX1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UHG3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AEL RAPPDKIAII GAGIGGTSAA YYLRQKFGKD VKIDLFEREE VGGRLATMMV QGQEYEAGGS VIHPLNLHMK RFVKDLGLSA VQASGGLLGI YNGETLVFEE SNWFIINVIK LVWRYGFQSL RMHMMWVEDVL DKFMRIYRYQ SHDYAFSSVE KLLHALGGDD FLGMLNRTLL ETLQKAGFSE KFLNEMIAPV MRVNYGQSTD INAFVGAVSL SCSDSGLWAV EGGNKLVCSSG LLQASKSNLI SGSVMYIEEK TKTKYTGNTPT KMYEVVYQIG TETRSDFYDI VLVATPLNRK MSNITFLNFD PPIEEFHQYY QHIVTTLVKG ELNTSIFSSR PIDKFGLNTV LTTDNSDLFI NSIGIVPSVR EKEDPEPSTD GTYVWKIFSQ ETLTKAQILK LFLSYDYAVK KPWLAYPHYK PPEKCPSSIL HDRLYYLNIGI ECAASAMEMS AIAAHNAALL AYHRWNGHTD MIDQDGLYEK LKTEL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PCYOX1
<b>Protein Names</b>	Recommended name: Prenylcysteine oxidase 1 EC= 1.8.3.5 Alternative name(s): Prenylcysteine lyase
<b>Expression Region</b>	28-505
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.