



Recombinant Human Zinc finger protein 559 (ZNF559)

Product Code	CSB-EP026839HU
Abbreviation	ZNF559
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BR84
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVAGWLTNYS QDSVTFEDVA VDFTQEEWTL LDQTQRNLYR DVMLENYKNL VAVDWESHIN TKWSAPQQNF LQGKTSSVVE MERNHFGEEL FDFNQCEKAL SEHSCLKTHR RTYFRKKTCE CNQCEKAFRK PSIFTLHKKT DIGEELPNCN QCETAFSQHL HLVCKKTSQN LHLVCKKHTT QEKPYKCSDC EKGLPSSSHL RECVRIYGGG RPYTHKEYVE TFSHSTALFV HMQTQDGEKF YECKACGKPF TESSYLQHL RTHSRVLPPIE HKKFGKAFAP SPDLAKHIRL RTRGKHYYVCN ECGKEFTCFS KLNIHIRVHT GEKPYECNKC GKAFDSSGL IKHRRHTTGE KPYECKECGK AFANSSHLTV HMRTHTGEKP YQCKEKGKAF INSSSFKSHM QTHPGVKPYD CQQCGKAFIR SSFLIRHLRS HSAERPFECE ECGKAFRYSS HLSQHKRIHT GERPYKCQKC GQAFSISGL TVHMRHTTGE RPFECQECGK AFTRSTYLIR HLRSHSVEKP YKECGQTFNS SSCLTECV
Source	E.coli
Target Names	ZNF559
Protein Names	Recommended name: Zinc finger protein 559
Expression Region	1-538
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.