



# Recombinant Human Zinc finger protein 559 (ZNF559)

<b>Product Code</b>	CSB-MP026839HU
<b>Abbreviation</b>	ZNF559
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BR84
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVAGWLTNYS QDSVTFEDVA VDFTQEEWTL LDQTQRNLYR DVMLENYKNL VAVDWESHIN TKWSAPQQNF LQGKTSSVVE MERNHFGEEL FDFNQCEKAL SEHSCLKTHR RTYFRKKTCE CNQCEKAFRK PSIFTLHKKT DIGEELPNCN QCETAFSQHL HLVCKKTSQN LHLVCKKHTT QEKPYKCSDC EKGLPSSSHL RECVRIYGGG RPYTHKEYVE TFSHSTALFV HMQTQDGEKF YECKACGKPF TESSYLQHL RTHSRVLPPIE HKKFGKAF AF SPDLAKHIRL RTRGKHYVCN ECGKEFTCFS KLNHIRVHT GEKPYECNKC GKAFDSSGL IKHRRTHTGE KPYECKECGK AFANSSHLTV HMRTHTGEKP YQCKECKGAF INSSSFKSHM QTHPGVKPYD CQQCGKAFIR SSFLIRHLRS HSAERPFECE ECGKAFRYSS HLSQHKRIHT GERPYKCQKC GQAFSISGL TVHMRTHTGE RPFECQECGK AFTRSTYLIR HLRSHSVEKP YKECGQTFNS SSCLTECV
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF559
<b>Protein Names</b>	Recommended name: Zinc finger protein 559
<b>Expression Region</b>	1-538
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.