



Recombinant Human Retbindin (RTBDN)

Product Code	CSB-BP020560HU
Abbreviation	RTBDN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BSG5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SRPLQARSQQ HHGLAADLGK GKLHLGAPCC PSEMDTTETS GPGNHPERCG VPSPECESFL EHLQRALRSR FRLRLLGVRQ AQLPCEELCQ AWFANCEDDI TCGPTWLPLS EKRGCPSCL TYGQTFADGT DLCRSALGHA LPVAAPGARH CFNISISAVP RPRPGRRGRE APSRRRSRSPR TSILDAAGSG SGSGSGSGP
Source	Baculovirus
Target Names	RTBDN
Protein Names	Recommended name: Retbindin
Expression Region	31-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.