



Recombinant Clupea harengus Cathepsin D (ctsd)

Product Code	CSB-BP006187CMS
Abbreviation	ctsd
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9DEX3
Product Type	Recombinant Protein
Immunogen Species	Clupea harengus (Atlantic herring)
Purity	>85% (SDS-PAGE)
Sequence	APTPETLKN YMDAQYYGEI GLGTPVQMFT VVFDTGSSNL WLPSIHCSFT DIACLLHHKY NGAKSSTYVK NGTEFAIQYG SGSLSGYLSQ DSCTIGDIVV EKQLFGAEIK QPGVAFIAAK FDGILGMAYP RISVDGVPPV FDMMSQKKV EQNVFSFYLN RNPDTEPGGE LLLGGTDPKY YTGDFNYVPV TRQAYWQIHM DGMSIGSQLT LCKDGCEAIV DTGTSLITGP PAEVRALQKA IGAIPLIQGE YMIDCKKVPT LPTISFNVGG KTYSLTGEQY VLKESQGGKT ICLSGLMGLE IPPPAGPLWI LGDVFQYY TVFDRESNRV GFAKST
Source	Baculovirus
Target Names	ctsd
Protein Names	Recommended name: Cathepsin D EC= 3.4.23.5
Expression Region	62-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>