



Recombinant Human Nuclear RNA export factor 2 (NXF2)

Product Code	CSB-EP016219HU-B
Abbreviation	NXF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZY0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MCSTLKKCGT YRTEVAECHD HGSTFQGRKK GGSSFRDNFD KRSCHYEHGG YERPPSHCQE NDGSVEMRDV HKDQQLRHTP YSIRCERRMK WHSEDEIRIT TWRNRKPPER KMSQNTQDGY TRNWFKVTIP YGIKYDKAWL MNSIQSHCSD RFTPVDHFHYV RNRACFFVQD ASAASALKDV SYKIYDDENQ KICIFVNHST APYSVKNKLLK PGQMEMLKLK MNKRYNVSQQ ALDLQNLRFDPDLMGRDIDI ILNRRNCMAA TLKIIERNFP ELLSLNLCNN KLYQLDGLSD ITEKAPKVKT LNLSKNKLES AWELGKVKGL KLEELWLEGN PLCSTFSDQS AYVSAIRDCF PKLLRLDGRE LSAPVIVDID SSETMKPCKE NFTGSETLKH LVLQFLQYY SIYDSGDRQG LLGAYHDEAC FSLAIPFDPK DSAPSSLCKY FEDSRNMKTL KDPYLKCELL RRTKRDIVDS LSALPKTQHD LSSILVDVWC QTERMLCFVS NGVFKEVEGQ SQGSVLAFTTRFIATPGSSS SLCIVNDELFRDASPQETQ SAFSIPVSTL SSSSEPSLSQ EQQEMVQAFS AQSGMKLEWS QKCLQDNEWN YTRAGQAFTM LQTEGKIPAE AFKQIS
Source	E.coli
Target Names	NXF2
Protein Names	Recommended name: Nuclear RNA export factor 2 Alternative name(s): Cancer/testis antigen 39 Short name= CT39 TAP-like protein 2 Short name= TAPL-2
Expression Region	1-626
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.