



# Recombinant Human Nuclear RNA export factor 2 (NXF2)

<b>Product Code</b>	CSB-MP016219HU
<b>Abbreviation</b>	NXF2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9GZY0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MCSTLKKCGT YRTEVAECHD HGSTFQGRKK GGSSFRDNFD KRSCHYEHGG YERPPSHCQE NDGSVEMRDV HKDQQLRHTP YSIRCERRMK WHSEDEIRIT TWRNRKPPER KMSQNTQDGY TRNWFKVTIP YGIKYDKAWL MNSIQSHCSD RFTPVDHFHYV RNRACFFVQD ASAASALKDV SYKIYDDENQ KICIFVNHST APYSVKNKLNK PGQMEMLKLTK MNKRYNVSQQ ALDLQNLRFDPDL MGRDIDI ILNRRNCMAA TLKIIERNFP ELLSLNLCNN KLYQLDGLSD ITEKAPKVKT LNLSKNKLES AWELGKVKGL KLEELWLEGN PLCSTFSDQS AYVSAIRDCF PKLLRLDGRE LSAPVIVDID SSETMKPCKE NFTGSETLKH LVLQFLQYY SIYDSGDRQG LLGAYHDEAC FSLAIPFDPK DSAPSSLCKY FEDSRNMKTL KDPYLKCELL RRTKRDIVDS LSALPKTQHD LSSILVDVWC QTERMLCFSV NGVFKEVEGQ SQGSVLAFTTR TFIATPGSSS SLCIVNDELFRD ASPQETQ SAFSIPVSTL SSSSEPSLSQ EQQEMVQAFS AQSGMKLEWS QKCLQDNEWN YTRAGQAFTM LQTEGKIPAE AFKQIS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	NXF2
<b>Protein Names</b>	Recommended name: Nuclear RNA export factor 2 Alternative name(s): Cancer/testis antigen 39 Short name= CT39 TAP-like protein 2 Short name= TAPL-2
<b>Expression Region</b>	1-626
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.