



# Recombinant Human Otoraplin (OTOR)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP017277HU  |
| <b>Abbreviation</b>      | OTOR  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9NRC9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | VHG IFMDRLASKK LCADDECVYT ISLASAQEDY NAPDCRFINV<br>KKGQQIYVYS KLVKENGAGE FWAGSVYGDG QDEMGVVGYP<br>PRNLVKEQRV YQEATKEVPT TDIDFFCE  |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | OTOR  |
| <b>Protein Names</b>     | Recommended name: Otoraplin Alternative name(s): Fibrocyte-derived protein<br>Melanoma inhibitory activity-like protein   |
| <b>Expression Region</b> | 18-128  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Target Details</b>    | This protein is secreted via the Golgi apparatus and may function in cartilage development and maintenance. A frequent polymorphism in the translation start codon of this gene can abolish translation and may be associated with forms of deafness. This gene is a member of the melanoma-inhibiting activity gene family. In addition, alternate polyA sites exist for this gene.                        |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |