



Recombinant Human Integrator complex subunit 9 (INTS9)

| | |
|--------------------------|--|
| Product Code | CSB-BP011767HU |
| Abbreviation | INTS9 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9NV88 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MKLYCLSGHP TLPCNVLKFK STTIMLDCGL DMTSTLNFLP LPLVQSPRLS NLPGWSLKDG NAFLDKELKE CSGHVFVDSV PEFCLPETEL IDLSTVDVIL ISNYHCMMAL PYITEHTGFT GTVYATEPTV QIGRLLMEEL VNFIERVPA QSASLWKNKD IQRLLPSPLK DAVEVSTWRR CYTMQEVNSA LSKIQLVGYS QKIELFGAVQ VTPLSSGYAL GSSNWIQSH YEKVSYVSGS SLLTTHPQPM DQASLKNSDV LVLTLGLTQIP TANPDGMVGE FCSNLALTVR NGGNVLVPCY PSGVIYDLLE CLYQYIDSAG LSSVPLYFIS PVANSSLEFS QIFAEWLCHN KQSKVYLPEP PFPHAELIQT NKLKHYPYSIH GDFSNDFRQP CVVFTGHPSL RFGDVVHFME LWGKSSLNTV IFTEPDFSYL EALAPYQPLA MKCIYCPIDT RLNFIQVSKL LKEVQPLHVV CPEQYQPPP AQSHRMDLMI DCQPPAMSYR RAEVLALPFK RRYEKIEIMP ELADSLVPME IKPGISLATV SAVLHTKDNK HLLQPPRPA QPTSGKKRKR VSDDVPDCKV LKPLLSGSIP VEQFVQTLEK HGFSDIKVED TAKGHIVLLQ EAETLIQIEE DSTHIICDND EMLRVRLRDL VLKFLQKF |
| Source | Baculovirus |
| Target Names | INTS9 |
| Protein Names | Recommended name: Integrator complex subunit 9 Short name= Int9 Alternative name(s): Protein related to CPSF subunits of 74 kDa Short name= RC-74 |
| Expression Region | 1-658 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.