



Recombinant Human Integrator complex subunit 9 (INTS9)

Product Code	CSB-EP011767HU-B
Abbreviation	INTS9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NV88
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKLYCLSGHP TLPCNVLKFK STTIMLDCGL DMTSTLNFLP LPLVQSPRLS NLPGWSLKDG NAFLDKELKE CSGHVFVDSV PEFCLPETEL IDLSTVDVIL ISNYHCMMAL PYITEHTGFT GTVYATEPTV QIGRLLMEEL VNFIERVPA QSASLWKNKD IQRLLPSPLK DAVEVSTWRR CYTMQEVNSA LSKIQLVGYS QKIELFGAVQ VTPLSSGYAL GSSNWIIQSH YEKVSYVSGS SLLTTHPQPM DQASLKNSDV LVLTLTQIP TANPDGMVGE FCSNLALTVR NGGNVLVPCY PSGVIYDLLE CLYQYIDSAG LSSVPLYFIS PVANSSLEFS QIFAEWLCHN KQSKVYLPEP PFPHAELIQT NKLKHYPYSIH GDFSNDFRQP CVVFTGHPSL RFGDVVHFME LWGKSSLNTV IFTEPDFSYL EALAPYQPLA MKCIYCPIDT RLNFIQVSKL LKEVQPLHVV CPEQYQPPP AQSHRMDLMI DCQPPAMSYR RAEVLALPFK RRYEKIEIMP ELADSLVPME IKPGISLATV SAVLHTKDNK HLLQPPRPA QPTSGKKRKR VSDDVPDCKV LKPLLSGSIP VEQFVQTLEK HGFSDIKVED TAKGHIVLLQ EAETLIQIEE DSTHIICDND EMLRVRLRDL VLKFLQKF
Source	E.coli
Target Names	INTS9
Protein Names	Recommended name: Integrator complex subunit 9 Short name= Int9 Alternative name(s): Protein related to CPSF subunits of 74 kDa Short name= RC-74
Expression Region	1-658
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.