



# Recombinant Human Spermine oxidase (SMOX)

<b>Product Code</b>	CSB-EP021844HU
<b>Abbreviation</b>	SMOX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NWM0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQSCESSGDS ADDPLSRGLR RRGQPRVVVI GAGLAGLAAA KALLEQGFTD VTVLEASSHI GGRVQSVKLG HATFELGATW IHGSHGNPIY HLAEANGLLE ETTDGERSVG RISLYSKNGV ACYLTNHGRR IPKDVVEEFS DLYNEVYNLT QEFFRHDKPV NAESQNSVGV FTREEVRNRI RNDPDDPEAT KRLKLAMIQQ YLKVESCESS SHSMDEVSLS AFGEWTEIPG AHHIIPSGFM RVVELLAEGI PAHVIQLGKP VRCIHWQAS ARPRGPEIEP RGEDHNDHT GEGGQGGEEP RGGRWDEDEQ WSVVVECEDC ELIPADHVIV TVSLGVLKRQ YTSFFRPGLP TEKVAAIHRL GIGTTDKIFL EFEEPFWGPE CNSLQFVWED EAESHTLTYP PELWYRKICG FVDLYPPER Y GHVLSGWICG EEALVMEKCD DEAVAEICTE MLRQFTGNPN IPKPRILRS AWGSNPYFRG SYSYTQVGSS GADVEKLAKP LPYTESSKTA PMQVLFSGEA THRKYYSTTH GALLSGQREA ARLIEMYRDL FQQGT
<b>Source</b>	E.coli
<b>Target Names</b>	SMOX
<b>Protein Names</b>	Recommended name: Spermine oxidase EC= 1.5.3.16 Alternative name(s): Polyamine oxidase 1 Short name= PAO-1 Short name= PAOh1
<b>Expression Region</b>	1-555
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The product of this gene is the polyamine oxidase. This enzyme potentially represents a new class of catabolic enzymes in the mammalian polyamine metabolic pathway capable of the efficient oxidation of polyamines. More than five transcript variants encoding four active isoenzymes have been identified for this gene, however, not all variants have been fully described. The characterized isoenzymes have distinctive biochemical characteristics and substrate specificities, suggesting the existence of additional levels of complexity in polyamine catabolism.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.