



Recombinant Human Zinc finger protein 586 (ZNF586)

Product Code	CSB-BP026869HU
Abbreviation	ZNF586
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NXT0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAAAALRAP AQSSVTFEDV AVNFSLEEWS LLNEAQRCLY RDVMLETLTL ISSLGCWHGG EDEAAPSKQS TCIHIYKDQG GHSGERPYPEC GEYRKLFFKNK SCLTEPRRDH KHRNVRTGER PYECSKYGKL FHQKPTLHIH ERFHTGQKTY ECSECGKSFH QSSSLLQRQT LHTRERPYPEC IECGKAFAEK SSLINHRKVH SGAKRYECNE CGKSFAYTSS LIKHRRHTG ERPYECSECG RSFAENSSLI KHLRVHTGER PYECVECGKS FRRSSSLLQH QRVHTRERPYPY ECSECGKSFS LRSNLIHHQR VHTGERHECG QCGKSF SRKS SLIIHLRVHT GERPYECSDC GKSFAENSSL IKHLRVHTGE RPYECIDCGK SFRHSSSFRR HQRVHTGMRP YK
Source	Baculovirus
Target Names	ZNF586
Protein Names	Recommended name: Zinc finger protein 586
Expression Region	1-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.