



# Recombinant Human Zinc finger protein 224 (ZNF224)

<b>Product Code</b>	CSB-MP026592HU
<b>Abbreviation</b>	ZNF224
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9NZL3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTTFKEAMTF KDVAVFTEE ELGLLDLAQR KLYRDVMLEN FRNLLSVGHQ AFHRDTFHFL REEKIWMMKT AIQREGNSGD KIQTEMETVS EAGTHQEWSF QQIWEKIASD LTRSQDLMIN SSQFSKEGDF PCQTEAGLSV IHTRQKSSQG NGYKPSFSDV SHFDFHQQLH SGEKSHTCDE CGKNFCYISA LRIHQRVHMG EKCYKCDVCG KEFSQSSHLQ THQRVHTGEEK PFKCVECCKG FSRRSALNVH HKLHTGEKPY NCEECGKAFI HDSQLQEHQR IHTGEEKPKFC DICGKSFCGR SRLNRHSMVH TAEKPFRCDT CDKSFRQRSA LNSHRMIHTG EKPYKCEECG KGFICRRDLY THHMVHTGEEK PYNCKECKGKS FRWASCLLKH QRVHSGEKPF KCEECGKGFY TNSQCYSHQR SHSGEKPYKC VECGKGYKRR LDLDFHQRVH TGEKLYNCKE CGKSFSRAPC LLKHERLHSG EKPFQCEECG KRFTQNSHLH SHQRVHTGEEK PYKCEKCGKG YNSKFNLDMH QKVHTGERPY NCKECKGSFG WASCLLKHQR LHSGEKPKFC EECGKRFTQN SQLHSHQRVH TGEKPYKCDE CGKGFWSST RLTHQRRHSR ETPLKCEQHG KNIVQNSFSK VQEKVHSVEK PYKCEDCGKG YNRRLNLDLMH QRVHMGEKTW KCRECDMCF S QASSLRLHQN VHVGEKP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF224
<b>Protein Names</b>	Recommended name: Zinc finger protein 224 Alternative name(s): Bone marrow zinc finger 2 Short name= BMZF-2 Zinc finger protein 233 Zinc finger protein 255 Zinc finger protein 27 Zinc finger protein KOX22
<b>Expression Region</b>	1-707
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.