



Recombinant Human Zinc finger protein 232 (ZNF232)

Product Code	CSB-EP026599HU-B
Abbreviation	ZNF232
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UNY5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAVSLTAAET LALQGTQGQE KMMMMGPKKE EQSCEYETRL PGNHSTSQEI FRQRFRHLRY QETPGPREAL SQLRVLCCEW LRPEKHTKEQ ILEFLVLEQF LTIPEELQS WVRGHHPKSG EEA VTVLEDL EGGLEPEPQV PGPAHGPAQE EPWEKKESLG AAQEALSIQL QPKETQPFPPK SEQVYLHFLS VVTEDGPEPK DKGSLPQPPI TEVESQVFSE KLATDTSTFE ATSEGTLELQ QRNPKAERLR WSPAQEEFR QMVVIHKEIP TGKKDHECSE CGKTFIYN SH LVVHQRVHSG EKPYKCSDCG KTFKQSSNLG QHQRIHTGEK PFECNECGKA FRWGAHLVQH QRIHSGEKPY ECNECGKAFS QSSYLSQHRR IHSGEKPFIC KECGKAYGWC SELIRHRRVH ARKEPSH
Source	E.coli
Target Names	ZNF232
Protein Names	Recommended name: Zinc finger protein 232 Alternative name(s): Zinc finger and SCAN domain-containing protein 11
Expression Region	1-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.