



Recombinant Human Zinc finger protein 148 (ZNF148)

Product Code	CSB-BP026544HU
Abbreviation	ZNF148
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UQR1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MNIDDKLEGL FLKCGGIDEM QSSRTMVVMG GVSGQSTVSG ELQDSVLQDR SMPHQEILAA DEVLQESEMR QQDMISHDEL MVHEETVKND EEQMETHERL PQGLQYALNV PISVKQEITF TDVSEQLMRD KKQIREPVDL QKKKKRKQRS PAKILTINED GSLGLKTPKS HVCEHCNAAF RTNYHLQRHV FIHTGEKPFQ CSQCDMRFIQ KYLLQRHEKI HTGEKPFRC D ECGMRFIQKY HMERHKRTHS GEKPYQCEYC LQYFSRTDRV LKHKRMCHEN HDKKNLNRCAI KGGLLTSEED SGFSTSPKDN SLPKKKRQKT EKKSSGMDKE SALDKSDLKK DKNDYLPYLS SSTKVKDEYM VAEYAVEMPH SSVGGSHLED ASGEIHPPKL VLKKINSKRS LKQPLEQNQT ISPLSTYEEES KVS KYAFELV DKQALLDSEG NADIDQVDNL QEGPSKPVHS STNYDDAMQF LKKKRYLQAA SNN SREYALN VGTIASQPSV TQAAVASVID ESTTASILES QALNVEIKSN HDKNVIPDEV LQTL LDHYSH KANGQHEISF SVADTEVTSS ISINSSEVPE VTPSENVGSS SQASSSDKAN MLQEYSKFLQ QALDRTSQND AYLNSPSLNF VTDNQTL PNQ PAFSSIDKQV YATMPINSFR SGMNSPLRTT PDKSHFGLIV GDSQHSFPFS GDET N HASAT STQDFLDQVT SQKKA EAQPV HQAYQMSSFE QPFRAPYHGS RAGIATQFST ANGQVNL RGP GTSAEFSEFP LVNVNDNRAG MTSSPDATTG QTFG
Source	Baculovirus
Target Names	ZNF148
Protein Names	Recommended name: Zinc finger protein 148 Alternative name(s): Transcription factor ZBP-89 Zinc finger DNA-binding protein 89
Expression Region	1-794
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.