



Recombinant Human ERBB receptor feedback inhibitor 1 (ERRFI1)

Product Code	CSB-EP007804HU
Abbreviation	ERRFI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UJM3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SIAGVAAQE IRVPLKTGFL HNGRAMGNMR KTYWSSRSEF KNNFLNIDPI TMAYSLNSSA QERLIPLGHA SKSAPMNGHC FAENGPSQKS SLPPLLIPPS ENLGPHEEDQ VVCGFKLTV NGVCASTPPL TPIKNPSLFCAPLCERGS RPLPPLPISE ALSLDDTDCE VEFLTSSDTD FLLEDSTLSD FKYDVPGRRS FRGCGQINYA YFDTPAVSAA DLSYVSDQNG GVPDPNPPPP QTHRRLRRSH SGPAGSFNKP AIRISNCCIH RASPNSDEDK PEVPPRVPIP PRPVKPDYRR WSAEVTSSTY SDEDRPPKVP PREPLSPSNS RTPSPKSLPS YLNGVMPPTQ SFAPDPKYVS SKALQRQNSE GSASKVPCIL PIIENGKKVS STHYYLLPER PPYLDKYEKF FREAETNGG AQIQPLPADC GISSATEKPD SKTKMDLGGH VKRKHLSYVV SP
Source	E.coli
Target Names	ERRFI1
Protein Names	Recommended name: ERBB receptor feedback inhibitor 1 Alternative name(s): Mitogen-inducible gene 6 protein Short name= MIG-6
Expression Region	2-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.