



Recombinant Human REST corepressor 1 (RCOR1), Partial

Product Code	CSB-EP019515HU
Abbreviation	RCOR1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UKL0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVEKGPEVSG KRRGRNNAAS SASAAAASAA ASAACASPAA TAASGAAASS ASAAAASAAA APNNGQNKSL AAAAPNGNSS SNSWEEGSSG SSSDEEHGGG GMRVGPQYQA VVPDFDPAKL ARRSQERDNL GMLVWSPNQN LSEAKLDEYI AIAKEKHGYN MEQALGMLFW HKHNIKSLA DLPNFTPFPD EWTVEDKVLV EQAFSFGKT FHRIQQMLPD KSIASLVKFY YSWKKTRTKT SVMDRHARKQ KRERESEDE LEEANGNNPI DIEVDQNKES KKEVPPTETV PQVKKEKHST QAKNRAKRKP PKGMFLSQED VEAVSANATA ATTVLRQLDM ELVSVKRQIQ NIKQTNALK EKLDGGIEPY RLPEVIQKCN ARWTTEEQLL AVQAIRKYGR DFQAISDVIG NKSVVQVKNF FVNYRRRFNI DEVLQEWAE HGKEETNGPS NQKPKVSPDN SIKMPEEEDE APVLDVRYAS AS
Source	E.coli
Target Names	RCOR1
Protein Names	Recommended name: REST corepressor 1 Alternative name(s): Protein CoREST
Expression Region	4-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.