



# Recombinant Human REST corepressor 1 (RCOR1), Partial

<b>Product Code</b>	CSB-YP019515HU
<b>Abbreviation</b>	RCOR1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UKL0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MVEKGPEVSG KRRGRNNAAS SASAAAASAA ASAACASPAA TAASGAAASS ASAAAASAAA APNNGQNKSL AAAAPNGNSS SNSWEEGSSG SSSDEEHGGG GMRVGPQYQA VVPDFDPAKL ARRSQERDNL GMLVWSPNQN LSEAKLDEYI AIAKEKHGYN MEQALGMLFW HKHNIKSLA DLPNFTPFPD EWTVEDKVLV EQAFSFGKT FHRIQQMLPD KSIASLVKFY YSWKKTRTKT SVMDRHARKQ KRERESEDE LEEANGNNPI DIEVDQNKES KKEVPPTETV PQVKKEKHST QAKNRAKRKP PKGMFLSQED VEAVSANATA ATTVLRQLDM ELVSVKRQIQ NIKQTNALK EKLDGGIEPY RLPEVIQKCN ARWTTEEQLL AVQAIRKYGR DFQAISDVIG NKSVVQVKNF FVNYRRRFNI DEVLQEWAE HGKEETNGPS NQKPKVSPDN SIKMPEEEDE APVLDVRYAS AS
<b>Source</b>	Yeast
<b>Target Names</b>	RCOR1
<b>Protein Names</b>	Recommended name: REST corepressor 1 Alternative name(s): Protein CoREST
<b>Expression Region</b>	4-485
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.