



# Recombinant Human Zinc finger protein 643 (ZNF643)

<b>Product Code</b>	CSB-YP026913HU
<b>Abbreviation</b>	ZNF643
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UJL9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLQQLLITLP TEASTWVKLR HPKAATERVA LWEDVTKMFK AEALLSQDAD ETQGESLESR VTLGSLTAES QELLTFKDVS VDFTQEEWGQ LAPAHRNLYR EVMLENYGNL VSVGCLSKP GVISQLEKGE EPWLMERDIS GVPSSDLKSK TKTKESALQN DISWEELHCG LMMERFTKGS SMYSTLGRIS KCNKLESQQE NQRMGKGQIP LMCKKTFTQE RGQESNRFEK RINVKSEVMP GPIGLPRKRD RKYDTPGKRS RYNIDLNVHS RSYTKMKTFE CNICEKIFKQ LIHLTEHMRI HTGEKPFRCCK ECGKAFSQSS SLIPHQRIHT GEKPYECKEC GKTRHPSSL TQHVRIHTGE KPYECRVCEK AFSQSIGLIQ HLRTHVREKP FTCKDCGKAF FQIRHLRQHE IIHTGVKPYI CNVCSKTFESH STYLTQHQRT HTGERPYKCK ECGKAFSQRI HLSIHQRVHT GVKPYECSHC GKAFRHDSSF AKHQRIHTGE KPYDCNECGK AFSCSSSLIR HCKTHLRNTF SNVV
<b>Source</b>	Yeast
<b>Target Names</b>	ZFP69B
<b>Protein Names</b>	Recommended name: Zinc finger protein 643
<b>Expression Region</b>	1-534
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.