



Recombinant Schizosaccharomyces pombe Uroporphyrinogen decarboxylase (hem12)

Product Code	CSB-MP025678SXV
Abbreviation	hem12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9USJ5
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKNDLILRAA KGEEVERPPV WIMRQAGRYL PEYHKLRKQ SFFEMCQTPE TACELTLQPV TRFKGLLDAA IIFSDILVIP QALGMQVVML EQKGPFFPKP LVVPEDIDLL EKTPNISAKL GYVMDAISLT REKLDGQVPL MGFSGAPWTI MAYMIEGGGS KTFAKAKSWL FRYPEASHKL LKIITDATVS YLIQQVYAGA QLLQIFDSWA GELSPEDFTE YAYPYLVRIC QEVKQHLKKK KRDEVPMIVF AKGAWYAIQ LCDSGYDVIG LDWTVSPKEA VRIRGNRRVT FQGNLDPNIL YGTREIIEAR TKEMIQDFGG GKQGYIINLG HGITPGVNPD DVRFFLEKCH QYGSA
Source	Mammalian cell
Target Names	hem12
Protein Names	Recommended name: Uroporphyrinogen decarboxylase Short name= UPD Short name= URO-D EC= 4.1.1.37
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.