



Recombinant Human Zinc finger protein 221 (ZNF221)

Product Code	CSB-EP026589HU-B
Abbreviation	ZNF221
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UK13
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MISPSLELLH SGLCKFPEVE GKMTTFKEAV TFKDVAVVFT EEELGLLDPA QRKLYRDVML ENFRNLLSVG NQPFHQDTFH FLGKEKFWKM KTTSQREGNS GGKIQIEMET VPEAGPHEEW SCQQIWEQIA SDLTRSQNSI RNSSQFFKEG DVPCQIEARL SISHVQQKPY RCNECKQSFS DVSVFDLHQQ SHSGEKSHTC GECGKSFCYS PALHIHQRVH MGEKCYKCDV CGKEFNQSSH LQTHQRVHTG EKPFKCGQCG KGFHSRSALN VHCKLHTGEK PYNCEECGKA FIHDSQLQEH QRIHTGEKPF KCDICGKSFR VRSRLNRHSM VHTGEKAFRC DTCGKNFRQR SALNSHSMVH IEEKPYKCEQ CGKGFICRRD FCKHQMVHTG EKPYNCKECG KTFRWSSCLL NHQQVHSGQK SFKCEECGKG FYTNSRRSSH QRSHNGEKPY NCEECGKDYK RRLDLEFHQR VHTGERPYNC KECGKSFGWA SCLLKHQRLH SGEKPFKCEE CGKRFTQSSQ LSHQTCHTG EKLYKCEQCE KGYNSKFNLD MHQRVHGGER PYNCKECGKS FGWASCLLKH QRLHSGEKPL KSGVWEEIYS EFTASFTSVS LCGRKA I
Source	E.coli
Target Names	ZNF221
Protein Names	Recommended name: Zinc finger protein 221
Expression Region	1-617
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.