



Recombinant Human Zinc finger protein 215 (ZNF215)

Product Code	CSB-EP026585HU-B
Abbreviation	ZNF215
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UL58
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MQPLSKLMAI SKPRNLSLRE QREVLRADMS WQQETNPVVE THDSEASRQK FRHFQYLKVS GPHEALSQW ELCLQWLRPE IHTKKQIIEI LVLEQFLAIL PEEVRTWVNL QHPNNSKDMV TLIEDVIEML EDEMPCKDS ALQMGSIKEK MKAGSRTGKP QEPVTFKDVV VEFKSKEEWGQ LDSAVKNLYR NVMLENFRNL NSLRKAHLLS KPFESLKLES KKKRWIMEKE IPRKTIFDMK SISGEESHG VIMTRLTESG HPSSDAWKGE NWLYRNQKKW DINLPQEAFI PETIYTEEED FECSENKKSF DINSVSSICA IQVGIPSRKG SPKCDKFKTY FKNLDSVKG QHSEYEGND LSLSTDIRHQ KSHTTMNSYE CYQCGKAFCR SSSLIRHQII HTGEKPYKCS ECGRFFNRRT NLTKHQKLHA EAKACTSNKC GKAFSKSEDS NNPTLHFGNN FYQCVNCGKS FNRSSSLIRH QMIHTGEKPF KCKECSKAFN RSSNLVKHQK LHTRDKS
Source	E.coli
Target Names	ZNF215
Protein Names	Recommended name: Zinc finger protein 215 Alternative name(s): BWSCR2-associated zinc finger protein 2 Short name= BAZ-2 Zinc finger protein with KRAB and SCAN domains 11
Expression Region	1-517
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.