



# Recombinant Human Zinc finger protein 215 (ZNF215)

<b>Product Code</b>	CSB-YP026585HU
<b>Abbreviation</b>	ZNF215
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UL58
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MQPLSKLMAI SKPRNLSLRE QREVLRADMS WQQETNPVVE THDSEASRQK FRHFQYLKVS GPHEALSQW ELCLQWLRPE IHTKKQIIEI LVLEQFLAIL PEEVRTWVNL QHPNNSKDMV TLIEDVIEML EDEMPCKDS ALQMGSIKEK MKAGSRTGKP QEPVTFKDVV VEFSKEEWGQ LDSAVKNLYR NVMLENFRNL NSLRKAHLLS KPFESLKLES KKKRWIMEKE IPRKTIFDMK SISGEESHG VIMTRLTESG HPSSDAWKGE NWLYRNQKKW DINLPQEAFI PETIYTEED FECSENKKSF DINSVSSICA IQVGIPSRKG SPKCDKFKTY FKNLDSVKG QHSEYEGND LSLSTDIRHQ KSHTTMNSYE CYQCGKAFCR SSSLIRHQII HTGEKPYKCS ECGRFFNRRT NLTKHQKLHA EAKACTSNKC GKAFSKSEDS NNPTLHFGNN FYQCVNCGKS FNRSSSLIRH QMIHTGEKPF KCKECSKAFN RSSNLVKHQK LHTRDKS</p>
<b>Source</b>	Yeast
<b>Target Names</b>	ZNF215
<b>Protein Names</b>	Recommended name: Zinc finger protein 215 Alternative name(s): BWSCR2-associated zinc finger protein 2 Short name= BAZ-2 Zinc finger protein with KRAB and SCAN domains 11
<b>Expression Region</b>	1-517
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.