



Recombinant Human Exocyst complex component 7 (EXOC7)

| | |
|--------------------------|--|
| Product Code | CSB-BP007886HU |
| Abbreviation | EXOC7 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9UPT5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MIPPQEASAR RREIEDKCLKQ EEETLSFIRD SLEKSDQLTK NMVSILSSFE SRLMKLENSI IPVHKQTENL QRLQENVEKT LSCLDHVISY YHVASDTEKI IREGPTGRLE EYLGSMKIQ KAVEYFQDNS PDSPELNKVK LLFERGKEAL ESEFRSLMTR HSKVVSPVLI LDLISGDDDL EAQEDVTLEH LPESVLQDVI RISRWLVEYG RNQDFMNVYY QIRSSQLDRS IKGLKEHFHK SSSSSGVPYS PAIPNKRKDT PTKKPVKRPG TIRKAQNLLK QYSQHGLDGK KGGSNLIPLE GLLPCTPRGG LPGPWINAAC VCAADISPGH EHDFRVKHLS EALNDKHGPL AGRDDMLDVE TDAYIHCVSA FVKLAQSEYQ LLADIPEHH QKKTFDLSLIQ DALDGLMLEG ENIVSAARKA IVRHDFSTVL TVFPILRHLK QTKPEFDQVL QGTAASTKNK LPGLITSMET IGAKALEDFA DNIKNDPDK EYNMPKDGTVH ELTSNAILFL QQLLDFQETA GAMLASQETS SSATSYSSEF SKRLLSTYIC KVLGNLQLNL LSKSKVYEDP ALSAIFLHNN YNYILKSLEK SELIQLVAVT QKTAERSYRE HIEQQIQTYQ RSWLKVTDYI AEKNLPVFQP GVKLRDKERQ I IKERFKGFN DGLEELCKIQ KAWAIPDTEQ RDRIRQAQKT IVKETYGAFL QKFGSVPFTK NPEKYIKYGV EQVGD MIDRL FDTSA |
| Source | Baculovirus |
| Target Names | EXOC7 |
| Protein Names | Recommended name: Exocyst complex component 7 Alternative name(s): Exocyst complex component Exo70 |
| Expression Region | 1-735 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.