



Recombinant Human Exocyst complex component 7 (EXOC7)

Product Code	CSB-EP007886HU-B
Abbreviation	EXOC7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UPT5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MIPPQEASAR RREIEDKCLKQ EEETLSFIRD SLEKSDQLTK NMVSILSSFE SRLMKLENSI IPVHKQTENL QRLQENVEKT LSCLDHVISY YHVASDTEKI IREGPTGRLE EYLGSMKIQ KAVEYFQDNS PDSPELNKVK LLFERGKEAL ESEFRSLMTR HSKVVSPVLI LDLISGDDDL EAQEDVTLEH LPESVLQDVI RISRWLVEYG RNQDFMNVYY QIRSSQLDRS IKGLKEHFHK SSSSSGVPYS PAIPNKRKDT PTKKPVKRPG TIRKAQNLLK QYSQHGLDGK KGGSNLIPLE GLLPCTPRGG LPGPWINAAC VCAADISPGH EHDFRVKHL S EALNDKHGPL AGRDDMLDVE TDAYIHCVSA FVKLAQSEYQ LLADIPEHH QKKTFDLSLIQ DALDGLMLEG ENIVSAARKA IVRHDFSTVL TVFPILRHLK QTKPEFDQVL QGTAASTKNK LPGLITSMET IGAKALEDFA DNIKNDPDK E YNMPKDGTVH ELTSNAILFL QQLLDFQETA GAMLASQETS SSATSYSSEF SKRLLSTYIC KVLGNLQLNL LSKSKVYEDP ALSAIFLHNN YNYILKSLEK SELIQLVAVT QKTAERSYRE HIEQQIQTYQ RSWLKVTDYI AEKNLPVFQP GVKLRDKERQ I IKERFKGFN DGLEELCKIQ KAWAIPDTEQ RDRIRQAQKT IVKETYGAFL QKFGSVPFTK NPEKYIKYGV EQVGDMIDRL FDTSA
Source	E.coli
Target Names	EXOC7
Protein Names	Recommended name: Exocyst complex component 7 Alternative name(s): Exocyst complex component Exo70
Expression Region	1-735
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.