



Recombinant Human Zinc finger protein 223 (ZNF223)

Product Code	CSB-BP026591HU
Abbreviation	ZNF223
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UK11
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTMSKEAVTF KDVAVVFTEE ELGLLDLAQR KLYRDVMLEN FRNLLSVGHQ PFHRDTFHFL REEKFWMMDI ATQREGNSGG KIQPEMKTFF EAGPHEGWSC QQIWEEIASD LTRPQDSTIK SSQFFEQGDA HSQVEEGLSI MHTGQKPSNC GKCKQSFSDM SIFDLPQQIR SAEKSHSCDE CGKSFCYISA LHIHQRVHLG EKLKCDVCG KEFSQSLHLQ THQRVHTGEK PFKCEQCGRG FRCRSALTVH CKLHMGEKHY NCEACGRAFI HDFQLQKHQR IHTGEKPFKC EICSVSFRLR SSLNRHCVVH TGKKPNSTGE YGKGFIRRLD LCKHQTIHTG EKPYNCKEKG KSFRRSSYLL IHQRVHTGEK PYKCDKCGKS YITKSGLDLH HRAHTGERPY NCDDCGKSFR QASSILNHKR LHCRKKPFKC EDCGKKLVYR SYRKDQQKNH SGENPSKCED CGKRYKRRNLN LDIIISLFLN DT
Source	Baculovirus
Target Names	ZNF223
Protein Names	Recommended name: Zinc finger protein 223
Expression Region	1-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.