



# Recombinant Human Zinc finger protein 223 (ZNF223)

<b>Product Code</b>	CSB-YP026591HU
<b>Abbreviation</b>	ZNF223
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UK11
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTMSKEAVTF KDVAVFTEE ELGLLDLAQR KLYRDVMLEN FRNLLSVGHQ PFHRDTFHFL REEKFWMMDI ATQREGNSGG KIQPEMKTFF EAGPHEGWSC QQIWEEIASD LTRPQDSTIK SSQFFEQGDA HSQVEEGLSI MHTGQKPSNC GKCKQSFSDM SIFDLPQQIR SAEKSHSCDE CGKSFCYISA LHIHQRVHLG EKLKCDVCG KEFSQSLHLQ THQRVHTGEK PFKCEQCGRG FRCRSALT VH CKLHMGEKHY NCEACGRAFI HDFQLQKHQR IHTGEKPFKC EICSVSFRLR SSLNRHCVVH TGKKPNSTGE YGKGFIRRLD LCKHQTIHTG EKPYNCKEKG KSFRRSSYLL IHQRVHTGEK PYKCDKCGKS YITKSGLDLH HRAHTGERPY NCDDCGKSFR QASSILNHR LHCRRKPFKC EDCGKKLVYR SYRKDQQKNH SGENPSKCED CGKRYKRRLN LDIILSLFLN DT
<b>Source</b>	Yeast
<b>Target Names</b>	ZNF223
<b>Protein Names</b>	Recommended name: Zinc finger protein 223
<b>Expression Region</b>	1-482
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.