



# Recombinant Human Zinc finger protein 225 (ZNF225)

<b>Product Code</b>	CSB-EP026593HU
<b>Abbreviation</b>	ZNF225
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9UK10
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTTLKEAVTF KDVAVVFTEE ELRLDLAQR KLYREVMLEN FRNLLSVGHQ SLHRDTFHFL KEEKFWMMET ATQREGNLGG KIQMEMETVS ESGTHEGLFS HQTWEQISSD LTRFQDSMVN SFQFSKQDDM PCQVDAGLSI IHVRQKPSEG RTCKKSFSDV SVLDLHQQLQ SREKSHTCDE CGKSFCYSSA LRIHQRVHMG EKLYNCDVCG KEFNQSSHLQ IHQRIHTGEK PFKCEQCGKG FSRRSGLYVH RKLHTGVKPH ICEKCGKAFI HDSQLQEHQR IHTGEKPFKC DICCKSFRSR ANLNRHSMVH MREKPFRCDT CGKSFGLKSA LNSHRMVHTG EKRYKCEECG KRFIYRQDLY KHQIDHTGEK PYNCKEKGKS FRWASGLSRH VRVHSGETTF KCEECGKGFY TNSQRYSHQR AHSGEKPYRC EECGKGYKRR LDLDFHQRVH RGEKPYNCKE CGKSFGWASC LLNHQRIHSG EKPFKCEECG KRFTQNSQLY THRRVHSGEK PFKCEECGKR FTQNSQLYSH RRVHTGVKPY KCEECGKGFN SKFNLDMHQR VHTGERPYNC KECGKSFSRA SSILNHNKRLH GDEKPFKCEE CGKRFTENSQ LHSHQRVHTG EKPYKCEKCG KSFRWASTHL THQRLHSREK LLQCEDCGKS IVHSSCLKDQ QRDQSGETKS KCEDCGKRYK RRLNLDTLLS LFLNDT
<b>Source</b>	E.coli
<b>Target Names</b>	ZNF225
<b>Protein Names</b>	Recommended name: Zinc finger protein 225
<b>Expression Region</b>	1-706
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.