



Recombinant Human Zinc finger protein 490 (ZNF490)

Product Code	CSB-MP026786HU
Abbreviation	ZNF490
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ULM2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRRNSSLSFQ MERPLEEQVQ SKWSSSQGRT GTGGSDVLQM QNSEHHGQSI KTQTDSISLE DVAVNFTLEE WALLDPGQRN IYRDVMRATF KNLACIGEKW KDQDIEDEHK NQGRNLRSPM VEALCENKED CPCGKSTSQI PDLNTNLETP TGLKPCDCSV CGEVFMHQVS LNRHMRSHT E QKPNECHEYG EKPHKCKEKG KTFTRSSSIR THERIHTG EK PYECKECKGKA FAF LFSFRNH IRIHTGETPY ECKECKGAFR YLTALRRHEK NHTG EKPYKC KQCGKAFIYY QPFLTHERTH TGEKPYECKQ CGKAFSCPTY LRSHEKHTTG EKPFVCRECG RAFFSHSSLR KHVKTHTGVQ PYTCKKCGEA FKSSSSCEVH ERTHFGEKPY ECKQCGKAFN SSSYLQLHER VHTG EKTYEC KECGKAFLYS THFRIHERTH TREKPYECKQ CGRVFIYFSH LRRHERSHTG VKPCECKQCG KAFTCLNSLK VHKRIHTGER PFQCRQCGKA FSYSKSLHVH ERTHSRQKP
Source	Mammalian cell
Target Names	ZNF490
Protein Names	Recommended name: Zinc finger protein 490
Expression Region	1-529
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.