



# Recombinant Mouse Mitogen-activated protein kinase 11 (Mapk11)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP013450MO-B  |
| <b>Abbreviation</b>      | Mapk11  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9WUI1  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Mus musculus (Mouse)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MSGPRAGFYR QELNKTWVEV PQLRQGLRPV GSGAYGSVCS<br>AYDARLRQKV AVKKLSRPFQ SLIHARRTYR ELRLLKHLKH ENVIGLLDVF<br>TPATSIEDFS EVYLVTTLMG ADLNNIVKCQ ALSDEHVQFL VYQLLRGLKY<br>IHSAGIIHRD LKPSNVAVNE DCELRILDFG LARQADEEMT GYVATRWYRA<br>PEIMLNWMHY NQTVDIWSVG CIMAELLQ GK ALFPGNDYID QLKRIMEVVG<br>TPSPEVLAKI SSEHARTYIQ SLPPMPQKDL SSVFHGANPL AIDLLGRMLV<br>LDSDQRVSAA EALAHAYFSQ YHDPDDEPEA EPYDESVEAK ERTLEEWKEL<br>TYQEVLSFKP LEPSQLPGTH EIEQ  |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | Mapk11  |
| <b>Protein Names</b>     | Recommended name: Mitogen-activated protein kinase 11 Short name= MAP kinase 11 Short name= MAPK 11 EC= 2.7.11.24 Alternative name(s): Mitogen-activated protein kinase p38 beta Short name= MAP kinase p38 beta Short  |
| <b>Expression Region</b> | 1-364   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Target Details</b>    | This protein is a member of the MAP kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation, and development. This kinase is most closely related to p38 MAP kinase, both of which can be activated by proinflammatory cytokines and environmental stress. This kinase is activated through its phosphorylation by MAP kinase kinases (MKKs), preferably by MKK6. Transcription factor ATF2/CREB2 has been shown to be a substrate of this kinase. |

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.