



Recombinant Nuclear RNA export factor 2 (nxf-2)

Product Code	CSB-MP016219CXY
Abbreviation	nxf-2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9XVS8
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MRGQNRRGYR NIEGRLSLSS HSSHSSPRQT HVTNLRPAEI KVVQDVVDGC YVATDDVLNL SNFSKNTEFV ERDMLMCLTK TRVMSVVLQH IGYKYPRISG ISFSNNRLCH LDHLSSLSSI SKFLKFLDLS HNQVLKVFVF LNFITFFSKI SSGEELKCLG TIPVETVFFE GNPVCEKFVQ CAEYANFIQK TFPKCSNLDG MEVEPKPDHN RIEQIIPFRN GYYGSDEVRT LVEEFIITYY KIYDGADGQQ TRKQLLDAYD TNNSTFTHTV VCLWDPIKFV MYPDSESYRM YLRTSHNVLN QEYFAANRAS RISHGAMDIV VALSRLPATI HLMDTFVVDV FLVSATLLGF TLHGTFRDGP SAIKPENTEE HDNYFTRTFM VAPRGEGKVA IVSDQLFISS MSKRREGDQYR MLVETATDID Q
Source	Mammalian cell
Target Names	nxf-2
Protein Names	Recommended name: Nuclear RNA export factor 2
Expression Region	1-421
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>