



Recombinant Human Zinc finger protein 835 (ZNF835)

Product Code	CSB-EP027077HU-B
Abbreviation	ZNF835
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2P0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEGLLSVALQ GAELEGNWKH EGQVEDLQEN QESCPEPEAV ACKGDPAGDS MQRERDEF SRI PRTISSPAAT QASVPDDSSS RRCSAPGES P KERHPDSRQR ERGGGPKKPW KCGDCGKA FS YCSAFILHQR IHTGEKPFAC PECGKA FSQS VHLTLHQ RTH TGEKPYACHE CGKA FSQGSY LASHWRHTG EKPHRCADCG KAFTRVHTLT QHRRVHTGER PYACAQCAKA FRNRSSLIEH QRIHTGEKPY ECSACA KAFR FSSALIRHQR IHTEEKPYRC GQCAKAF AQI AHLTQHRRVH TGEKPYTCQD CGALFSQSAS LAEHRRIHTG EKP YACGQCA KAFTQVSHLT QHQ RTH TGER PYPCHDCGKR FSNRSHLLQH RLVHTGERPY RCLQCGA AFS HVSSLIEHQK IHTGERPYKC GECGKA FSQGSY SSLALHQ RTH TGERPYTCPE CGKA FSNRSY LIQHIVHTG EKP YECGSGCG KA FSFSSALI RHQRTHADSS GRLCPAPTPD STPGLSQGGE TCQQGCPGRN PRGPAED
Source	E.coli
Target Names	ZNF835
Protein Names	Recommended name: Zinc finger protein 835
Expression Region	1-537
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.