



# Recombinant Human Zinc finger protein 256 (ZNF256)

<b>Product Code</b>	CSB-EP026615HU-B
<b>Abbreviation</b>	ZNF256
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9Y2P7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAAAELTAPA QGIVTFEDVA VYFSWKEWGL LDEAQKCLYH DVMLENLTLT TSLGGSGAGD EEAPYQQSTS PQRVSQVRIP KALPSPQKTN PCEICGPVLR QILHLVEHQG THHGQKLYTD GACRKQLQFT AYLHQHQKQH VGQKHFRSNG GRDMFLSSCT FEVSGKPFTC KEVGKDFLVR SRFLQQAAH TRKKSNRTKS AVAFHSVKNH YNWGECVKAF SYKHVRVQHQ GDLIRERSYM CSECGKSFST SCSLSDHLRV HTSEKPYTCG ECGKSYRQSS SLITHRRIHT GVRPHQCDEC GKLFNRKYDL LIHQRVHTGE RPYKCSECGK SFSHSSLIT HQRHTGMRP YECSECGKSF IHSSSLITHQ RVHTGTRPYM CSECGKSFSQ SCHLIKHRRL HIGEGPYECS ECGKLFTYRS RFFQHQRVHT GVRSHCEHEC GKLF SRKFDL IVHERVHTGE RPYECSECGK SFTCKSYLIS HWKVHTGARP YECGECGKSF THSSTLLQHQ RVHTGERPYE CNECGKFFSQ SSSLIRHRRS HTGERPYECS ECWKSFSNHS SLVKHRRVHT GERPYECSEC GKSF SQSSNL TNHQRHSGE RPYECSDCGK FFTFNSNLLK HQNVHKG
<b>Source</b>	E.coli
<b>Target Names</b>	ZNF256
<b>Protein Names</b>	Recommended name: Zinc finger protein 256 Alternative name(s): Bone marrow zinc finger 3 Short name= BMZF-3
<b>Expression Region</b>	1-627
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.