



Recombinant Rat Runt-related transcription factor 2 (Runx2)

Product Code	CSB-MP020594RA
Abbreviation	Runx2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Z2J9
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MLHSPHKQPQ NHKCGANFLQ EDSKEALVFK WLISAGHYQP PRPTESVSAL STVHAVIFKA ASSIYNRGHK FYLEKKGGM ASNSLFSAVT PCQQSFFWNK TLPVAFKVVA LGEVPDGTVV TVMAGNDENY SAELRNASAV MKNQVARFND LRFVGRSGRG KSFTLTITVF TNPPQVATYH RAIKVTVDGP REPRRHRQKL DDSKPSLFSR RLSDLGRI
Source	Mammalian cell
Target Names	Runx2
Protein Names	Recommended name: Runt-related transcription factor 2 Alternative name(s): Core-binding factor subunit alpha-1 Short name= CBF-alpha-1 Osteoblast-specific transcription factor 2 Short name= OSF-2
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is a member of the RUNX family of transcription factors and encodes a nuclear protein with an Runt DNA-binding domain. This protein is essential for osteoblastic differentiation and skeletal morphogenesis and acts as a scaffold for nucleic acids and regulatory factors involved in skeletal gene expression. The protein can bind DNA both as a monomer or, with more affinity, as a subunit of a heterodimeric complex. Mutations in this gene have been associated with the bone development disorder cleidocranial dysplasia (CCD). Transcript variants that encode different protein isoforms result from the use of alternate promoters as well as alternate splicing.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.