



Recombinant Mouse Eukaryotic translation initiation factor 2-alpha kinase 1 (Eif2ak1)

Product Code	CSB-YP007510MO
Abbreviation	Eif2ak1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Z2R9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MLGGSSVDGE RDTDDDAAGA VAAPPAIDFP AEVSDPKYDE SDVPAELQVL KEPLQQPTFP FLVANQLLLV SLLEHLSHVH EPNPLHSKQV FKLLCQTFIK MGLSSFTCS DEFSSLRLHH NRAITHLMRS AKERVRQDPC QDNSYMQKIR SREIAFEAQT SRYLNEFEEL AILGKGGYGR VYKVRNKLDG QHYAIKKILI KSATKTDCMK VLREVKVLG LQHPNIVGYH TAWIEHVHVV QPQDRVPIQL PSLEVLSEQE GDRDQGGVKD NESSSSIVFA ELTPEKEKPF GESEVKNENN NLVSYTANLV VRNSSESESS IELQEDGLTD LSVRPVVRHQ LPLGHSSELE GNFTSTDESS EGNLNLGQT EVRYHMLHI QMQLCELSLW DWITERNKRS REYVDEAACP YVMASVATKI FQELVEGVFY IHNMGIVHRD LKPRNIFLHG PDQQVKIGDF GLACADIQN ADWTNRNGKG TRTHTSRVGT CLYASPEQLE GSQYDAKSDM YSLGVILLEL FQPFGTEMER ATVLTGVRTG RIPESLSKRC PVQAKYIQLL TGRNVSQRPS ALQLLQSELF QTTGNVNLT L QMKIIEQEKE IEELKKQLSL LSQDRGLKR
Source	Yeast
Target Names	Eif2ak1
Protein Names	Recommended name: Eukaryotic translation initiation factor 2-alpha kinase 1 EC= 2.7.11.1 Alternative name(s): Heme-controlled repressor Short name= HCR Heme-regulated eukaryotic initiation factor eIF-2-alpha kinase Heme-regul
Expression Region	1-619
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein acts at the level of translation initiation to downregulate protein synthesis in response to stress. The encoded protein is a kinase that can be inactivated by hemin. Two transcript variants encoding different isoforms have been found for this gene.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.