



# Recombinant Human Zinc finger protein 844 (ZNF844)

<b>Product Code</b>	CSB-EP027086HU-B
<b>Abbreviation</b>	ZNF844
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q08AG5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDLVAFEDVA VNFTQEEWSL LDPSQKNLYR EVMQETLRNL ASIGEKWKDQ NIEDQYKNPR NNLRSLLGER VDENTEENHC GETSSQIPDD TLNKKTSPGV KCESSVCGE VFGHSSLNR HIRADTAHKP SEYQEYQGEP YKCQQRKKAF RCHPSFQMQE KAHTGEKLYD CKECGKTFIS HSSIQRHMIM HNGDGTYKCK FCGKACPCLS IYLIHERVHT GEKPYKCKQC GKAFSYSTSL QIHERTHTGE KPYECKECKG AFGSPNSLYE HRRTHTGEKP YECKQCGKAF RWFHSFQIHE RTHSEEKAYE CTKCGKAFKC PSYLCRHEVT HSGKKPCECK QCGKALSYN FQRHMKMTR MRPYKCKTVE KPLILPVRFE DMKELTLERN LMNASTVVKP SIVVPFTIM KGLTLERNPM NVSSVVKPSF LPLPFDIMKG LTLERNRMSV SNVGKPSDLP HTFKCMEGLT LKRNP MNVSS VVKPSFFPLP FDIMKGLTLE RNPMSVSNVG KPSHLPHTFK CMKGLTLESN CMNLNNVKKP LDLSETFKFM KRHTLERNPI RNMEKHSTIS LPFKYMQQCT EDRMP MNVKS VTKHSYLP RS FEYMQEHTLE RNP MNVRNAE KRSIIFLLCV YTKGCTLERN HINVRIVGKH SVCLVPFVDI KGLTLE
<b>Source</b>	E.coli
<b>Target Names</b>	ZNF844
<b>Protein Names</b>	Recommended name: Zinc finger protein 844
<b>Expression Region</b>	1-666
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.