



Recombinant Escherichia coli O6:K15:H31 Choline dehydrogenase (betA)

Product Code	CSB-EP005335EGY
Abbreviation	betA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TKW1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQFDYIIIGA GSAGNVLATR LTEDPNTTVL LLEAGGPDYR FDFRTQMPAA LAFPLQGKRY NWAYETEPEP FMNNRRMECG RGKGLGGSSL INGMCYIRGN ALDLNWAQE PGLNWSYLD CLPYRKAET RDVGENDYHG GDGPVSVITS KPGVNPLFEA MIEAGMQAGY PRTDDLNGYQ QEGFGPMDRT VTPHGRRAST SRGYLDQAKS RPNLTIRTHA MTDHIIFDGK RAVGVEWLEG DSTIPTRAAA NKEVLLCAGA IASPQILQRS GVGNAELLAE FDIPLVHELP GVGENLQDHL EMYLQYECKE PVSYPALQW WNQPRIGAEW LFGGTGVGAS NHFEAGGFIR SREEFAWPNI QYHFLPVAIN YNGSNAVKEH GFQCHVGSMR SPSRGHVRIK SRDPHQHPGI LFNYSHEQD WQEFRDAIRI TREIMHQPAL DQYRGREISP GVECQTDEQL DEFVRNHAET AFHPCGTCKM GYDEMAVVDG EGRVHGLEGL RVVDASIMPQ IITGNLNATT IMIGEKIADM IRGKEALPRS TAGYFVANGM TVRAKK</p>
Source	E.coli
Target Names	betA
Protein Names	Recommended name: Choline dehydrogenase Short name= CDH Short name= CHD EC= 1.1.99.1
Expression Region	1-556
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.