



Recombinant Escherichia coli O6:K15:H31 Ornithine carbamoyltransferase (argI)

Product Code	CSB-BP017270EGY
Abbreviation	argI
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0T9E0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MATSLKNRNF LKLLDYTPAE IQYLIDLAIN LKAAKKSNE KQTLVGKNIA LIFEKSSTRT RCAFEEVAADF QGAQVTYIGP SGSQIGHKES MKDTARVLGR MYDGEYRGY GQNIVEELGE FAGVPVWNGL TNEFHPTQIL ADLMTMLEHA PGKTLPELSF AYLGDARNNM GNSLMVGAAG MGMDIRLVAP KSFWPDEALV TQCREIASVT GARITLTEDV EEGVYDVDFL YTDVWVSMGE PKEAWAERVS LMAPYQINQQ VITATRNPEV KFMHCLPAFH NEHTTVGREI EMAYGSKDWK LLKRFLNRPT LLFSMKQKTG CTPLKR
Source	Baculovirus
Target Names	argI
Protein Names	Recommended name: Ornithine carbamoyltransferase Short name= OTCase EC= 2.1.3.3
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.