



Recombinant Human Zinc finger protein 169 (ZNF169)

Product Code	CSB-BP026553HU
Abbreviation	ZNF169
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14929
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSPGLLTTRK EALMAFRDVA VAFTQKEWKL LSSAQRTRYR EVMLENYSHL VSLGIAFSKP KLIEQLEQGD EPWRENEHL LDLCPEPRTE FQPSFPHLVA FSSSQLLRQY ALSGHPTQIF PSSSAGGDFQ LEAPRCSSEK GESGETEGPD SSLRKRPSRI SRTFFSPHQG DPVEWVEGNR EGGTDLRLAQ RMSGGSDTM LKGADTSESG AVIRGNYRLG LSKKSSLSFSH QKHHVCPECG RGFCQRSOLI KHQRTHTGEK PYLCPECGRR FSQKASLSIH QRKHSGEKPY VCRECGRHFR YTSSLTNHNR IHSGERPFVC QECGRGFRQK IALLHQRTTH LEEKPFVCPE CGRGFCQKAS LLQHQSSTHG ERPFLCLECG RSFRQQSLLL SHQVTHSGEK PYVCAECGHS FRQKVTLIRH QRTHTGEKPY LCPQCGRGFS QKVTLIGHQR THTGEKPYLC PDCGRGFGQK VTLIRHQRTTH TGEKPYLCPK CGRAFQFKSL LTRHQRTTHSE EELYVDRVCG QGLGQKSHLI SDQRTHSGEK PCICDECGRG FGFKSALIRH QRTHSGEKPY VCRECGRGFS QKSHLHRHRR TKSGHQLLPQ EVF
Source	Baculovirus
Target Names	ZNF169
Protein Names	Recommended name: Zinc finger protein 169
Expression Region	1-603
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.