



Recombinant Human Zinc finger protein 266 (ZNF266)

Product Code	CSB-BP026623HU
Abbreviation	ZNF266
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14584
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLENYKNLAT VGYQLFKPSL ISWLEQEESR TVQRGDFQAS EWKVQLKTKELALQQDVLGE PTSSGIQMIG SHNGGEVSDV KQCGDVSSEH SCLKTHVRTQ NSENTFECYL YGVDFLTLHK KTSTGEQRSV FSQCGKAFSL NPDVVCQRTC TGEKAFDCSD SGKSFINSH LQGHLRTHING ESLHEWKECG RGFHSTDLA VRIQTHRSEK PYKCKECKGK FRY SAYLNIH MGTHTGDNPY ECKECKGAFT RSCQLTQHRK THTGEKPYKC KDCGRAFTVS SCLSQHMKIH VGEKPYECKE CGIAFTRSSQ LTEHLKHTA KDPFECKICG KSFRNSSCLS DHFRIHTGIK PYKCKDCGKA FTQNSDLTKH ARTHSGERP Y ECKECKKAFA RSSRLSEHTR THTGEKPFEC VKCGKAFAIS SNLSGHLRIH TGEKPFECLE CGKAFTHSSS LNNHMRTHSA KKPFTCMECG KAFKFPTCVN LHMRIHTGEK PYKCKQCGKS FSYSNSFQLH ERTHTGEKPY ECKECKKAFA SSSSFRNHER RHADERLSA
Source	Baculovirus
Target Names	ZNF266
Protein Names	Recommended name: Zinc finger protein 266 Alternative name(s): Zinc finger protein HZF1
Expression Region	1-549
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.