



# Recombinant Human Zinc finger protein 266 (ZNF266)

<b>Product Code</b>	CSB-MP026623HU
<b>Abbreviation</b>	ZNF266
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q14584
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MLENYKNLAT VGYQLFKPSL ISWLEQEESR TVQRGDFQAS EWKVQLKTKELALQQDVLGE PTSSGIQMIG SHNGGEVSDV KQCGDVSSEH SCLKTHVRTQ NSENTFECYL YGVDFLTLHK KTSTGEQRSV FSQCGKAFSL NPDVVCQRTC TGEKAFDCSD SGKSFINSH LQGHLRTHING ESLHEWKECG RGFHSTDLA VRIQTHRSEK PYKCKECKGK FRY SAYLNIH MGTHTGDNPY ECKECKGAFT RSCQLTQHRK THTGEKPYKC KDCGRAFTVS SCLSQHMKIH VGEKPYECKE CGIAFTRSSQ LTHELKTHTA KDPFECKICG KSFRNSSCLS DHFRIHTGIK PYKCKDCGKA FTQNSDLTKH ARTHSGERP Y ECKECKKAFA RSSRLSEHTR THTGEKPFEC VKCGKAFAIS SNLSGHLRIH TGEKPFECLE CGKAFTHSSS LNNHMRTHSA KKPFTCMECG KAFKFTCVN LHMRIHTGEK PYKCKQCGKS FSYSNSFQLH ERTHTGEKPY ECKECKKAFA SSSSFRNHER RHADERLSA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF266
<b>Protein Names</b>	Recommended name: Zinc finger protein 266 Alternative name(s): Zinc finger protein HZF1
<b>Expression Region</b>	1-549
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.