



Recombinant Human Interferon-related developmental regulator 2 (IFRD2)

Product Code	CSB-MP011070HU
Abbreviation	IFRD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12894
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAWNPSRRP VWQGGAPRED GGARGVWLPS SGQVSAQRTG RRLVGGLEPTP TGSLTPRPPR PVPGMPRARK GNTLRKGGQR RGGGARSSAQ ADGSSDDEA ASEARSTASE CPSLLSTTAE DSLGGDVVDE QGQQEDLEEK LKEYVDCLTD KSAKTRQGAL ESLRLALASR LLPDFLLERR LTLADALEKC LKKGKGEEQA LAAAVLGLLC VQLGPGPKGE ELFHSLQPLL VSVLSDSTAS PAARLHCASA LGLGCYVAAA DIQDLVSCLA CLESVFSRFY GLGGSSTSPV VPASLHGLLS AALQAWALL TICPSTQISH ILDRQLPRLP QLLSESVNL RIAAGETIAL LFELARDLEE EFVYEDMEAL CSVLRTLATD SNKYRAKADR RRQRSTFRAV LHSVEGGECE EEIVRFGFEV LYMDSWARHR IYAAFKEVLG SGMHHHLQNN ELLRDIFGLG PVLLLDATAL KACKVPRFEK HLYNAAAFKA RTKARSVRD KRADIL
Source	Mammalian cell
Target Names	IFRD2
Protein Names	Recommended name: Interferon-related developmental regulator 2 Alternative name(s): Protein SKMC15
Expression Region	1-506
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.