



# Recombinant Human Splicing factor 1 (SF1)

<b>Product Code</b>	CSB-EP021120HU
<b>Abbreviation</b>	SF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q15637
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ATGANATPL DFPSKRRKRS RWNQDTMEQK TVIPGMPTVI PPGLTREQER AYIVQLQIED LTRKLRTGDL GIPPNPEDRS PSPEPIYNSE GKRLNTREFR TRKKLEEERH NLITEMVALN PDFKPPADYK PPATRVSDKV MIPQDEYPEI NFVGLLIGPR GNTLKNIEKE CNAKIMIRGK GSVKEGKVGR KDGQMLPGED EPLHALVTAN TMENVKKAVE QIRNILKQGI ETPEDQNDLR KMQLRELARL NGTLREDDNR ILRPWQSSET RSITNTTVCT KCGGAGHIAS DCKFQRPGDP QSAQDKARMD KEYLSLMAEL GEAPVPASVG STSGPATTPL ASAPRPAAPA NNPPPPSLMS TTQSRPPWMN SGPSESRPYH GMHGGGGPGGP GGGPHSFPHP LPSLTGGHGG HPMQHNPNGP PPPWMQPPPP PMNQGPHPPG HHGPPMDQY LGSTPVGSGV YRLHQGKGM PPPPMGMMPP PPPPPSGQPP PPPSGPLPPW QQQQQQPPPP PPSSSMASST TPLPWQQNTT TTTTSAGTGS IPPWQQQQA AAASPGAPQM QGNPTMVPLP PGVQPPLPPG APPPPPPPP GSAGMMYAPP PPPPPPMDPS NFVTMMGMGV AGMPPFGMPP APPPPPPQN
<b>Source</b>	E.coli
<b>Target Names</b>	SF1
<b>Protein Names</b>	Recommended name: Splicing factor 1 Alternative name(s): Mammalian branch point-binding protein Short name= BBP Short name= mBBP Transcription factor ZFM1 Zinc finger gene in MEN1 locus Zinc finger protein 162
<b>Expression Region</b>	2-639
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.